# **Assay Reproducibility**

## **Experimental variance of a BioFlux adhesion assay**

### **Experimental Design**

The microfluidic channels of a 24-well BioFlux 200 plate (Figure 1) were coated with human fibronectin (10  $\mu$ g/mL) for 1 hour. Channels were washed with RPMI + 10 mM HEPES for 10 minutes at 37°C. Cells were added to the inlet A wells, perfused for 5 minutes at 2 dyn/cm², and washed with PBS from inlet B wells for 10 minutes at 10 dyn/cm². Images were captured beginning at the outlet end of the viewing window.

Figure 1





Figure 1. BioFlux 24-well plate channels.

Microfluidic channels viewed from beneath a well plate. Inlet well are labeled "A" and "B" and outlet wells are labeled "O".

#### Figure 2



**Figure 2. Experimental Images.** Images from 3 experiments on the same plate were used for counting adherent cells. Images were taken at the same viewing location of each microfluidic channel.

#### **Results**

The average number of cells in all fields was 118 (STD=11). This represents a 9% variance across all populations captured and a 4% variance across the channels, which is in line with the normal variance of samples of cells in suspension.

	Channel 1	Channel 2	Channel 3
Avg. Cell Count	114	117	122
Standard deviation	10	13	12

For more information on BioFlux Technology, visit cellmicrosystems.com

