



# LAGO X

## Industry Leading Performance

Preclinical Optical Imaging | BLI, FLI & X-Ray

### Unrivalled Sensitivity

36% more sensitive when comparing minimum detectable radiance for BLI to leading competitor in its class. Patented LED excitation yields 100x more light on target for a stronger fluorescence signal with minimal background, resulting in superior quantitation.

### Ultimate in High Throughput


Industry leading 25cm x 25cm FoV allows for bioluminescence, fluorescence & X-ray imaging of up to 10 mice

### Robust & Reliable Workhorse

Solid state air-cooled cameras provide ultra cold -90°C absolute cooling without liquids or gas. No leaks & less pieces & parts to break.

### Free Analysis Software

Aura Software enables you to capture and analyze your data quickly, easily & accurately. **Easy mode** - just 3 clicks to acquisition at maximum sensitivity without compromising image quality. **Kinetics** - automatic, real-time acquisition of BLI kinetic curves.



To learn more about our family of instruments, visit [www.spectralinvo.com](http://www.spectralinvo.com)