

# Cytek<sup>®</sup> cFluor<sup>™</sup>


## TBMNK Kit, 8 Color

Turnkey Solution for Basic Immunoprofiling


### Cytek cFluor<sup>™</sup> RUO Immunoprofiling Kit

Designed by Cytek scientists to enable researchers to quickly perform basic immunophenotyping and analyze human B cell, T cell, monocytes, and NK cell populations, this ready-to-use cFluor TBMNK Kit is fully optimized for best performance on the Cytek Aurora and Northern Lights full spectrum systems.

The cFluor TBMNK Kit can also be used as a backbone to expand to a broader panel of your desire.



# 8-Color Sample Data



Violet		Blue		Red	
Specificity	Fluorochrome	Specificity	Fluorochrome	Specificity	Fluorochrome
CD3	cFluor V420	CD8	cFluor B515	CD16	cFluor R668
CD14	cFluor V450	CD19	cFluor BYG710	CD56	cFluor R720
CD45	cFluor V547			CD4	cFluor R780

Cytek dye cFluor™ R780 and its conjugates were developed, manufactured and commercialized by BioLegend, Inc. under the trademark APC/Fire™ 750. Fluor conjugated antibody manufactured and supplied by BioLegend, Inc.

cFluor™ V547, cFluor™ B515, cFluor™ R668 and cFluor™ R720 are equivalent to CF®405L, CF®488A, CF®647 and CF®700 respectively, manufactured and provided by Biotium, Inc. under an Agreement between Biotium and Cytek (LICENSEE). The manufacture, use, sale, offer for sale, or import of the product is covered by one or more of the patents or pending applications owned or licensed by Biotium. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel.

Find out more at [www.cytekbio.com/collection/reagents](http://www.cytekbio.com/collection/reagents)

Please contact [sales@cytekbio.com](mailto:sales@cytekbio.com) for any questions

For Research Use Only. Not for use in diagnostic procedures.

© 2021 Cytek Biosciences, Inc. All rights reserved. All trademarks are the properties of Cytek Biosciences Inc. and its subsidiaries unless otherwise specified.

N9-20074 Rev. A

