


Simoa Planar Array Human MMP Panel 1

Catalog Number 85-0007

MMP-3, MMP-9

The Simoa Planar Array Human Matrix Metalloproteinase (MMP) Panel 1 assay is a multiplex immunoassay designed for the Quanterix SP-X™ Imaging and Analysis system. The assay measures human MMP-3, MMP-9 in serum and plasma samples. Representative normal healthy sample readings are listed in the graph below along with the limit of detection (LoD) and limits of quantitation (LLoQ and ULoQ) for each biomarker in the assay. Simoa assay kits are available in 1-plex assay formats in addition to the multiplex format.



Human MMP Panel 1		
	MMP-3	MMP-9
LoD (pg/ml)	3.9	15.52
LLoQ (pg/ml)	9.76	48.8
ULoQ (pg/ml)	10,000	50,000